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Ver. 1.22

Table of Contents

1. LUMITESTER C-110 Utility Software.....	3
2. System requirements	3
3. Installation.....	4
4. Connecting to the LUMITESTER C-110	5
5. How to use	6
5-1 Transferring the measurement data from the LUMITESTER.....	6
(1) Transferring the data	6
(2) Saving the data	6
5-2 Creating a calibration curve.....	7
(1) Creating a calibration curve by measuring standard solutions	7
(2) Registering a calibration curve by direct input.....	11
6. Command description	12
6-1 File menu	12
6-2 Edit menu	12
6-3 Communication menu	12
6-4 Cal.Curve menu	12
6-5 Setup menu	13
7. Troubleshooting	13

Copyrights

Copyright © 2012 by Kikkoman Corporation. All Rights Reserved.

The contents of this manual are subject to change without prior notice.

Reproduction, transmission or copying of this manual, or saving for information retrieval, in part or in whole, as well as translation into other languages is prohibited without the permission of Kikkoman Corporation and Kikkoman Biochemifa Company.

Reproduction or use of this software for the purpose of printing, saving or retrieval, in part or in whole, is prohibited without the permission of Kikkoman Corporation.

Trademarks

- Windows is a registered trademark of Microsoft Corporation in the United States and/or other countries.
- Pentium and Celeron are trademarks or registered trademarks of Intel Corporation or its subsidiary in the United States and/or other countries.
- LabVIEW and NI-VISA are registered trademarks of National Instruments Corporation

Copyright (C) 2000 Microsoft corporation. All Rights Reserved.

Copyright (C) 2007 National Instruments corporation. All Rights Reserved.

1. LUMITESTER C-110 Utility Software

This software allows you to transfer the data acquired with the LUMITESTER C-110 to a PC (personal computer). You can then save the data in the PC or print it out as needed. You can also easily create a calibration curve using the Calibration Curve Creation Wizard and register it in the LUMITESTER.

2. System requirements

Hardware

- Intel Pentium/Celeron Family 600 MHz or faster processor
(1.5 GHz or faster processor recommended)
- 128 MB or more of RAM (512 MB or more of RAM recommended)
- 256 MB or more of available hard disk space
- DVD or CD-ROM drive

OS and corresponding version of the utility software

- Version 1.22 : Windows XP/Vista/7 (32bit/64bit)*

* Windows 7 compatibility

Our company performed verification of Windows 7 and our products of operation independently. Although we have taken all possible measures against verification of operation, a different result may occur depending on the environment of PC in use.

3. Installation

Start Windows and then insert the supplied CD into the CD drive of your PC.
The installer automatically starts. Follow the instructions that appear on the screen.


* In case of Windows Vista or 7:

When an auto-play file selection window is opened, please select and run a "C110-Utility_Installer.bat" file.

When a user account control window is opened, please choose "permit."

* In case the auto-play does not start:

Select the icon of the install CD in My Computer and right-click to open a context menu as shown in the right figure. Select and install utility software for the appropriate OS and language version.

- Uninstallation

When you uninstall this software, please perform uninstallation from the Add / Remove program of the control panel.

* The National Instruments folder created at the time of installation:

This software is using LabVIEW Runtime Engine and NI-VISA Runtime, which are products of National Instruments. These files are saved in the folder, "C:\Program Files\National Instruments". The main program of this application software is saved in a different folder.

- LabVIEW Runtime Engine

This is a software module, which is needed when performing application software.

- NI-VISA Runtime

NI-VISA Runtime is a module for using VISA I/O standard practice of NI-VISA, which National Instruments mounted.

* About NATIONAL INSTRUMENTS software license agreement

The "NATIONAL INSTRUMENTS SOFTWARE LICENSE AGREEMENT" displayed at the time of installation is a common license agreement for both of "Lumitester C-110 utility software" and a development environment "LabVIEW."

Moreover, "Software" described in the license agreement refers to the development environment "LabVIEW", and "Lumitester C-110 utility software" corresponds to the "Authorization Application".

4. Connecting to the LUMITESTER C-110

Check that the LUMITESTER C-110 and PC are off.

Connect the dedicated RS-232C cable (supplied) to the LUMITESTER C-110.

Connect the other end of the dedicated RS-232C cable to the serial connector (RS-232C) on your PC.

If your PC does not have a serial connector, then use a USB-serial conversion cable (USB-RS232C conversion cable) to connect to the USB port on your PC.

* Caution

If a cable is disconnected during communication between PC and the Lumitester, there is a possibility that an error may occur in the PC. The cable should not be disconnected until the communication is completed and the Lumitester is shut down.

Recommended cable

Make a choice between ① and ②

① USB-serial conversion adapter

Model : SRC06USM/SRC06USB Manufacturer : Buffalo Kokuyo Supply Inc.

RS232C straight cable (1.5m)

② Commercial USB-serial conversion adapter which is suitable for your OS.

*Caution!

Windows7 has two different OS (32 bit version / 64 bit version).

Please make sure your OS version before you choose an adapter.

5. How to use

5-1 Transferring the measurement data from the LUMITESTER

(1) Transferring the data

From the Communication menu, select "Transfer All Data" or "Transfer Every Time".


Transfer All Data

Selecting this menu command transfers the data stored in the LUMITESTER all at a time.

Transfer Every Time

Selecting this menu command transfers the result each time measurement is made with the LUMITESTER.

The measurement data is displayed as shown below after it is transferred.

	DATA#	RLU	Cal. Value	Time Stamp	Cal. Curve a	Cal. Curve b	COEF. (pg)	ATP/CELL
1	1	21,370		03/16/2010 13:49				
2	2	22,939	1.62E-08 M	03/16/2010 13:50	7.07E-13	-0.00E-00		
3	3	23,977	859.45 pg	03/16/2010 13:50	7.07E-13	-0.00E-00	5.07E+10	
4	4	22,429	1.59E+07 CELLS	03/16/2010 13:51	7.07E-13	-0.00E-00		1.00E-18
5	5	22,131	1.56E+07 CELLS	03/16/2010 13:51	7.07E-13	-0.00E-00		1.00E-18
6								
7								
8								
9								
10								
11								
12								
13								
14								
15								
16								
17								
18								
19								
20								

The RLU values, calibration values, calibration curve and coefficients used for calibration, etc. are displayed.

To stop the "Transfer Every Time" operation, click "Transfer Every Time" again on the Communication menu.

* Other menus cannot be selected during the "Transfer Every Time" operation.

* The number of data which can be displayed is 1000 data.

(2) Saving the data

From the File menu, select "Save As...", and name and save the data.

* The following two formats are used to save the data.

Dedicated file format

Data of this format can be displayed and printed only with this software.

Data is saved in this format when "Save As..." is selected.

CSV file format

This is a text format file with each value separated by a comma, and can be used with spreadsheet software.

Select "Export to CSV File" to save the data in CSV format.

5-2 Creating a calibration curve

There are two ways to create a calibration curve.

Measure Standard ATP

Selecting this menu command creates and registers a calibration curve by measuring standard solutions with known concentrations.

Input Directly

Selecting this menu command allows you to directly enter and register a calibration curve.

(1) Creating a calibration curve by measuring standard solutions

Form the Cal.Curve menu, select "Create" → "Measure Standard ATP".

The "Wizard 1/3 Input concentration of standard ATP solution" appears, as shown below.

Wizard to create the calibration curve. 1/3 Input concentration of standard ATP solution.

Input concentration of standard ATP solution and click the Next button.

Standard ATP	Concentration(M)
A	BLANK
B	<input type="text" value="1e-12"/>
C	<input type="text" value="3e-12"/>
D	<input type="text" value="6e-12"/>
E	<input type="text" value="1e-11"/>

1.00E-15 ~ 1.00E-4
Please input it by exponential value.
You will measure only the item which you input.

Cancel Next >>

Enter the concentrations of standard solution B to E.

(A is fixed to BLANK.)

To obtain a calibration curve, one or more standard solutions with known concentrations are required in addition to the BLANK solution (zero concentration).

After entering the conditions, click the "Next" button.

The "Wizard 2/3 Measurement of the standard solution" appears, as shown below.

Wizard to create the calibration curve. 2/3 Measurement of the standard solution...

Please measure standard solution with LUMITESTER.
Data is input to the position of the cursor.

Standard ATP	Concentration(M)	RLU
A	BLANK	0
B	1.00E-12	219
C	3.00E-12	589
D	6.00E-12	1,502
E	1.00E-11	2845

Cancel << Back Next >>

Specify the cell to make measurement.

Set a standard solution in the LUMITESTER and start measurement.


Measurement data is then displayed in the box where the cursor is located.

Make measurements of other standard solutions until data is entered in all empty boxes.

Then click the "Next" button.

* The measured values can also be directly input.

The "Wizard 3/3 Setting of the calibration curve" appears, as shown below.

A calibration curve calculated from the measurement data is displayed.

Click the "Registration" button to register the calibration curve in the LUMITESTER.

To display the calibration curve registered in the LUMITESTER, select "Show Cal.Curve" from the "Cal.Curve" menu.

The re-calculation of a calibration curve, and a graph display function

Re-calculation function

Put check mark in the "Fix b to 0" check box and click a re-calculation button, the calibration curve, which is fixed b to 0, is re-calculated.

R²

The coefficient of determination of an actual measurement and the prepared calibration curve is displayed. When the correlation is higher, this value becomes closer to 1.

Graph display

Click a graph button, and then a calibration curve will be displayed in the graph display window.

(1) Graph scale

The display scale of a graph can be changed. A linear scale or a log scale can be selected.

(2) Graph palette

Cursor button

It is not used.

Expansion / reduction button


Click the button will display the selection window of the range specification method. Choose the range specification method, drag or click a portion to be expanded in the graph, and perform scaling.

Scroll button

A graph can be scrolled, if a graph is dragged after clicking a button.

(2) Registering a calibration curve by direct input

From the Cal.Curve menu, select "Input Directly". The following dialog box appears.

Input a value of the calibration curve.

Input a value of the calibration curve.

model: $y = a \cdot x + b$

a	b(M)	[1E-17 ≤ a ≤ 1E-10]
7.07E-13	0.00E+0	[-1E-8 ≤ b ≤ 0]

The software sets the calibration curve to LUMITESTER.

OK Cancel

Specify the cell and enter a value.

After entering all values, click the "OK" button.

6. Command description

6-1 File menu

Command	Select to:
Open	Open a data file.
Close	Close the data file.
Save As...	Name and save the currently open data.
Export to CSV File	Save the data in CSV format which can be opened with spreadsheet software.
Print Preview	Display a print preview.
Print	Print the currently open data using a printer.
Exit	Quit this software.

6-2 Edit menu

Command	Select to:
Copy	Copy the specified data.
Select All	Select all data.
Clear All	Delete the displayed data.

6-3 Communication menu

Command	Select to:
Transfer All Data	Transfer the data stored in the LUMITESTER all at a time.
Transfer Every Time	Transfer the result each time measurement is made with the LUMITESTER.

6-4 Cal.Curve menu

Command	Select to:
"Create" → "Measure Standard ATP"	Create and register a calibration curve by measuring standard ATP solutions with known concentrations.
"Create" → "Input Directly"	Directly enter and register a calibration curve.
"Show Cal.Curve"	Display the calibration curve registered in the LUMITESTER.

6-5 Setup menu

Command	Select to:
Select Show Data	Set the data display items.
Date & Time Setup	Transfer the PC date/time to the LUMITESTER.
Show LUMITESTER Setup	Display or print the various settings for the LUMITESTER.

7. Troubleshooting

• Cannot communicate with the LUMITESTER C-110.

An error message "A communication error occurred" appears.

- i. When a communication error cannot be solved, please make re-connection in the following procedures.
 1. Turn off the Lumitester C-110.
 2. Close the utility software.
 3. Reconnect the dedicated RS-232C cable (supplied), which has connected PC and the Lumitester C-110, and turn on the Lumitester C-110.
 4. Start the utility software again.

If communication at the time of utility software starting is performed normally, the serial number of a Lumitester in use will be displayed on the title bar of the upper part of window.

- ii . Please communicate with Lumitester C-110 after measurement is completed.

• Cannot use commands.

"Transfer Every Time" on the Communication menu is being run.

Select "Transfer Every Time" again to terminate its operation.

When data is not input into a table, "Print Preview" and "Print" menu cannot be used.

Contact information

Kikkoman Biochemifa Company

2-1-1, Nishi-Shinbashi, Minato-ku, Tokyo 105-0003, Japan

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Copyright (C) 2012 by KIKKOMAN Corporation. All Rights Reserved.

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